

LewTM MDI 3.0 mm Multi-Platform Implant

o-ball



Platform Switching Collar

Provides a smooth interface for the soft tissue, creating a better emergence profile as well as hindrance for bacterial migration.

Labyrinth Groove Banding (Patent Pending)

Hinders cellular migration as well as assisting in cortical maintenance.

Reverse Buttress Threading

Allows ease of penetration into the bone and improves stability.

Cog LockTM (Patent Pending)

Anti-rotation barbs for added stability.



Multi-Platform

Used with Lew MDI O-Ball Housings and 3.0 mm Abutment Copings.

Progressive Compression Threading

Threads are gradually closer together along the length of the implant in order to provide a vertical compression effect and consequently added stability.

Optimized Tip for tactile feedback

Recommended Drilling Sequence for 3.0 mm implants

| Bone Type | Tissue Punch | Marking Drill | Drill 1.2 mm | Drill 1.5 mm | Drill 2.0 mm |
|-------------|--------------|---------------|--------------|-----------------|-----------------|
| Hard - I | 2.0 mm | • | | • Full Depth | • Half-Depth |
| Medium - II | 2.0 mm | • | • | | |
| Soft - III | 2.0 mm | • | • | | |

For Type I bone, a full depth 1.5 mm osteotomy should be drilled first. Then a 2.0 mm half-depth osteotomy should be drilled prior to implant insertion.

For Type II & III bone, a half-depth 1.2 mm osteotomy should be drilled.

As a general rule the depth should not exceed one-third to one-half the threaded length of the chosen MDI. The desired depth is the minimal amount that allows one to begin the auto-advancement of the implant into bone.



RESEARCH CENTER & LEGACY SALES
19 West 34th Street, Suite 916
New York, NY 10001

PARK PLAZA OFFICES
2401 N. Commerce, Suite E
Ardmore, OK 73401

212-736-3765


Toll Free: 800-243-7372

Facebook: www.facebook.com/ParkDental

www.parkdentalresearch.com

Instruments & Accessories

Site Preparation



MARK-1

Surgical Marker

TP-20

2.0 mm Tissue Punch

Surgical Drills & Drill Extender



MDT-11

Marking Drill

DRT-12

1.2 mm Surgical Drill

DRT-15

1.5 mm Surgical Drill

DRT-20

2.0 mm Surgical Drill

3.0 mm Multi-Platform Implant

Titanium



3.0 mm


| Lengths | Product Code |
|---------|--------------|
| 10 mm | 3010-CT |
| 13 mm | 3013-CT |
| 15 mm | 3015-CT |
| 17 mm | 3017-CT |

3.0 mm Head Analog



ACAN-3


Implant Drivers & Extender



CODRV-4

Torque Adapter/Extender


4 mm



DRVRT-204

O-Ball Mini Multi-Purpose Driver

2.04 mm



DRVRT-204L

O-Ball Mini Multi-Purpose Driver Long

2.04 mm

Wrenches



RATC-00

Park DentalTM Ratchet Wrench

RATC-80

Continuous Feedback Torque Wrench 15-80 Ncm

Includes CODRV-4

3.0 mm Impression Cap



OIC-3

3.0 mm Temporization Cap



OTC-3

3.0 mm Waxing Coping



OWC-3

Block Out Shim



BOS-1

O-Ball Housings



OMH-1



OMHP-1

O-Ball Housings



OMH-2



OMHP-2

O-Ring



OR-1

Micro O-Ring



OR-2